

Custom System Details:

Shell Type

Standard

Low Heel Cup
 Left Right
 Heel cup height is approx. 1/4 - 3/8

Medial Flange
 Left Right
 Heel Cup Height is approx. 1/2 - 5/8

Lateral Flange
 Left Right

Medial & Lateral
 Left Right

Roberts Plate
 Left Right

U.C.B.L.
 Left Right

Gait Plate
 Promote Out Toe
 Promote In Toe
 Left Right

Note: All height and widths are based on each individual foot. Approximations are based on the average foot.

Shell Dimensions

Width:

Narrow - Shell is contained between 1st and 5th met. heads
 Regular - Shell dichotomizes the 1st and 5th met. heads
 Wide - Shell encompasses the 1st and 5th met. heads

Standard

Arch Height:

To Cast - Determined during posting to balance the foot
 High - 1/8 of plaster fill
 Medium - 3/8 of plaster fill
 Low - 1/2 of plaster fill

Standard

Heel Cup Height:

Low - 3/8
 Medium - 5/8
 High - 3/4

Standard

Heel Raise

Left Right
 mm mm

Gait

Pronation Supination

Shell Padding

Padding may be added over the shell (under vinyl top only)

Standard

None
 1/16 Poron®
 1/8 Poron®

Posting Instructions (System I)

Post According To Positive Model (Standard)
 No Posting (neutral shell only)
 Forefoot Only (No measurements needed)
 Post According To Measurements

Left Varus Right Varus
 Left Varus Right Varus
 RearFoot

Left Varus Right Varus
 Left Varus Right Varus
 ForeFoot

Modifications To Shell

Cut Out 1st Ray In Shell Left Right
 Cut Out 5th Ray In Shell Left Right
 Reinforce Arch Left Right

Modifications To Metatarsal Bar (System II) Standard: 1/8 Thick 3° Varus Wedge

Changes By Degree:

Left Right
 Varus
 Valgus

Changes By Shape:

Left Right Left Right Left Right

Draw In Any Other Changes:

Left Right
 5 4 3 2 1 1 2 3 4 5
 Metatarsal Heads Metatarsal Heads

Modifications To Orthotic (System II)

Medial Bars Lateral Bars
 Left Right Left Right

X-13 Heel Lift

Male Female Shoe Size _____
 Height: _____ mm

Accommodation Chart:

Right Foot **Left Foot**

Mortons Extension

X for R or L Neuroma bar

X for R or L Accommodations

1/8" depth standard or

H Hole

Write in Desired Depth on selected Heads

X for Soft U-pad

small medium large
x for met. pad size and height

X for R or L Heel Spur pad

PLANTAR VIEW

Lateral Spur Pad Medial Spur Pad Center Spur Pad

Drill Out With Soft Plug